



Information below is for all technology levels of S Series products: 11, 9, 7, and 5.

Battery Life Expectancy Chart

Standard Products	Life Expectancy (in days)	Batteries Required (in 12 months)
-------------------	---------------------------	-----------------------------------

Size 13

S Series BTE	9 - 12	30 - 40
--------------	--------	---------

Size 312

S Series RIC 40	6 - 9	40 - 61
S Series RIC 50	6 - 9	40 - 61
S Series RIC AP 60	6 - 9	40 - 61
S Series RIC AP 71	4 - 7	52 - 92

Custom Products	Life Expectancy (in days)	Batteries Required (in 12 months)
-----------------	---------------------------	-----------------------------------

Size 13

S Series ITE	13 - 17	21 - 28
S Series HS/ITC	13 - 17	21 - 28
S Series CIC	13 - 17	21 - 28

Size 312

S Series ITE	7 - 10	36 - 52
S Series HS/ITC	7 - 10	36 - 52
S Series CIC	7 - 10	36 - 52

Size 10

S Series HS/ITC	5 - 7	52 - 73
S Series CIC	5 - 7	52 - 73

All estimates are based on a 16-hour day.